

Vaginitis Multi-Test Kit (Dry Chemoenzymatic Method)

Intended Use

This 6-parameter at-home vaginal health test gives you a simple way to track vaginal wellness markers like vaginal pH, healthy bacteria balance, signs of inflammation, and infections including yeast infections, BV, Aerobic Vaginitis (AV), and Trichomoniasis (Trich). Understanding your results helps you take action toward better health without frequent doctor visits. With easy instructions and a color chart, this test helps conveniently monitor your vaginal wellness.

Sample Requirements

Avoid sex, douching, vaginal creams, or suppositories for 24 hours before testing. Do not test if you are on your period, are using estrogen therapy, recently took antibiotics, or have misused hormones.

Testing Steps

Please follow the steps on the back of the included color card.

Interpretation of Results

pH

pH 3.8–4.4: Perfect. Your vagina's healthy "sour zone" (like mild vinegar). It's nature's way to block bad germs and protect good bacteria.

pH 4.6: (Mildly High) Common causes: sex fluids, period ending, normal body changes, or douching. Your body can often fix this!

Slight risk: Germs grow easier. Watch for: Odor, unusual discharge, itching.

pH 4.8: (Moderately High) Frequent cause: Bacterial Vaginosis (BV). Other: STIs (like Trich), douching, hormone changes.

pH 5.4 : (Much Higher) Typically indicates: Active BV or Trich. (Less common: Low hormones after menopause).

Hydrogen Peroxide (H₂O₂) = Levels of Good Bacteria

Negative ("–") Result: This is a GOOD result. It means you have plenty of good bacteria (lactobacilli) working to protect your vagina. These good bacteria naturally make hydrogen peroxide (H₂O₂), which acts like a cleaning agent to keep your vagina healthy and prevent infections.

Moderate ("±") Result: This means you have some good bacteria making small amounts of cleaning agents (H₂O₂), but not enough for full protection. Your vagina might need extra care to stay balanced and avoid infections.

Positive ("+") Result: means your vagina lacks enough good bacteria that naturally fight germs. Without enough of these protectors, your vagina becomes less protected. This makes you more likely to get infections like BV (bacterial vaginosis). Having plenty of these good bacteria lowers your BV risk.

Leukocyte (LE) = Inflammation (White Blood Cells)

("-" and "±") Result: YOU ARE GOOD!

Positive ("+", "++", "+++") Result means your vagina shows signs of inflammation — like a “warning light” telling you something might be wrong, but not what exactly. If you also test positive for specific infections like BV, trichomoniasis (“trich”), or yeast, this strongly confirms you need treatment.

If all **other tests are negative but LE is positive**, this could mean: a **very early infection** where germ levels are too low for other tests to detect yet (LE reacts first), non-infection irritation (from soaps, douching products, lubricants, condoms, hormonal changes like menopause or breastfeeding, or friction from sex or tampons), or other untreated issues like STIs (chlamydia/gonorrhea, not covered here) or general bacterial overgrowth (not BV/AV).

What to do next for your possible result? Track symptoms like unusual discharge, itching/burning, or pain/odor, avoid irritants like smelly soaps or douches, and consider regular at-home checks to see how your inflammation changes over time.

Sialidase (SNA) = BV

Negative ("–") Result: YOU ARE GOOD!

A “Possible BV” (“±”) Result: means your test was borderline—not clearly normal but not definitely BV. This can happen because: 1) You might be in the **early or healing stages of BV**, during which bacterial levels are changing; 2) Your vaginal bacteria are in a **mixed state** (some good, some bad); 3) **Recent activities** like sex, douching, or your menstrual cycle temporarily affected the test.

A Positive (“+”) Result strongly means you have **BV** (bacterial vaginosis), as the test detects specific chemicals produced by BV-causing bacteria like Gardnerella and Prevotella.

What to do next for your possible BV result?

Track symptoms:

- Watch for **fishy odor or thin gray discharge** — these suggest BV
- Note new discomfort during/after sex or periods

Optimize retest timing:

- **Avoid douching** — it harms your vagina’s natural protection
- **Wait a week** after antibiotics or your period ends
- **Retest** between periods for best results.
- **Avoid sex and lubricants** for 2 days before retesting

β -Glucuronidase (β -G) = Aerobic Vaginitis (AV)

Negative ("−") Result: YOU ARE GOOD!

A "**±**" (possible AV) Result means that the test detected borderline elevated levels of chemical signs produced by aerobic bacteria (AV) — above normal but not clearly positive. This outcome may occur during early infection, recovery, or due to sample collection issues, recent antibiotic use, or timing related to your period.

A **Positive ("+") Result** means **Aerobic Vaginitis (AV)**, an infection caused by harmful bacteria. While AV and Bacterial Vaginosis (BV) both cause itching/burning and may occur together, AV usually shows **yellow discharge with vaginal redness/swelling**, whereas BV features thin gray-white discharge with fishy odor but no redness/swelling. Importantly, some women—especially during pregnancy—have no AV symptoms; early testing is vital since untreated AV may cause early birth.

What to do next for your possible AV result?

- **Tracking symptoms daily**—note discharge (yellow) color, irritation, or swelling.
- **Timing retests carefully:** After your period ends or finishing antibiotics.
- **Temporarily avoiding irritants** like douches or plastic-feeling underwear.

N-acetylglucosaminidase (NAG) = Yeast or Trich

Negative ("−") Result: YOU ARE GOOD!

A "**Possible Trich or Yeast**" ("**±**") Result means your test shows borderline infection signs—not weak enough to be negative but not strong enough for a clear positive—with your vagina's sourness (acidity) in the middle zone. This could mean a very mild infection starting or healing. Sometimes blood, sex fluids, or recent douching affect the test; normal events like your period, recent sex, or using sprays/lubes can also change sourness (acidity) temporarily—watch for changes but don't worry yet.

(**NAG**) is a marker both trich and yeast share, so it can't tell them apart alone—but combined with your vaginal pH, it helps figure out whether you have trich or a yeast infection:

Positive ("+") and pH is high (4.8 or above), it's more likely to be **trich**. (Trich makes your pH higher.)

Positive ("+") and pH is low (4.6 or below), it's more likely to be a **yeast infection**.

What to do next for your "possible trich or yeast" result?

Retest in 5–7 days—avoid sex, douches, or vaginal creams for 24 hours beforehand, and wait until at least 3 days after your period ends.

Test Principle

Unlike UTI tests that only check for general problems, this test detects key vaginal health issues by identifying unique ‘chemical fingerprints’ (specific enzymes) in your vaginal fluid—precisely checking for healthy bacteria balance, Bacterial Vaginosis (BV), Aerobic Vaginitis (AV), yeast infections, Trichomoniasis (Trich), plus vaginal pH levels. It works by detecting enzymes only produced by these specific conditions, giving you more accurate results.

Storage & Expiration

Store at 36°F–86°F (2°C–30°C) in a dry place; use within 12 months. Once opened, test strips must be used within 1 hour. Always use before the expiration date shown on the package.

Testing Steps

Medical Interpretation Required: **This test does not replace professional medical diagnosis.**

Blood Specimens: Blood samples may yield inaccurate results. Use with extreme caution.

Handling Test Strips: Do not touch the reaction membrane (test area).

Single Use Only: For one-time use only. Discard after testing.

Testing Procedure: Ensure complete immersion of all reaction zones in sample fluid.

Safety Precautions: Keep out of children’s reach. Do not ingest.

Materials

Materials Provided

• Test Cassettes • Package Insert • Swab • Sample Tube • Color Chart

Materials Required But Not Provided

• Timer

Symbols

Item	Remark	Location	Item	Remark	Location
	LOT No.	Package		Manufacturer	User Manual Package
	Shelf Life	Package		Authorized Representative	User Manual Package
	Temperature: 2–30C	Package		CE Mark	User Manual Package
	Not Reuse	Package		Consult instructions for use	Package
	Date of Manufacture	Package		Avoid the Sun	Package
	In Vitro Diagnostic Medical Device	User Manual Package		Keep Dry	Package

